

| <b>SUND</b>                                                                                                                                             |     | <i>Rates as of 1/11-2025</i> |              |               |
|---------------------------------------------------------------------------------------------------------------------------------------------------------|-----|------------------------------|--------------|---------------|
| <b>Flow cytometer</b>                                                                                                                                   |     | <b>DIY/h</b>                 | <b>Min/h</b> | <b>Full/h</b> |
| Fortessa (3 laser)                                                                                                                                      | DKK | 304                          | 414          | 609           |
| CytoFLEX S (4 laser)                                                                                                                                    | DKK | 304                          | 414          | 609           |
| Fortessa (5 laser)                                                                                                                                      | DKK | 357                          | 467          | 661           |
| Fortessa X20 (5 laser)                                                                                                                                  | DKK | 357                          | 467          | 661           |
| FACSymphony A5                                                                                                                                          | DKK | 409                          | 519          | 714           |
| ImageStream                                                                                                                                             | DKK | 330                          | 441          | 635           |
| <b>Cell sorter</b>                                                                                                                                      |     |                              |              |               |
| FACS ARIA-II                                                                                                                                            | DKK | 430                          | 540          | 735           |
| FACS Melody                                                                                                                                             | DKK | 420                          | 530          | 724           |
| FACSymphony S6                                                                                                                                          | DKK | 483                          | 593          | 787           |
| <b>NGS instruments (booked time)</b>                                                                                                                    |     | <b>DIY/h</b>                 | <b>Min/h</b> | <b>Full/h</b> |
| 10X Chromium (run time app. 30 min)                                                                                                                     | DKK | 992                          | 1.102        | 1.296         |
| 10X Chromium X (run time app. 30 min)                                                                                                                   | DKK | 1.045                        | 1.150        | 1.335         |
| 10X Visium CytAssist (run time app. 30 min)                                                                                                             | DKK | 945                          | 1.050        | 1.235         |
| Rhapsody (run time app. 2-3 hours)                                                                                                                      | DKK | 278                          | 388          | 577           |
| Thermocycler                                                                                                                                            | DKK | 89                           | 199          | 393           |
| Tapestation (run time app. 30 min)                                                                                                                      | DKK | 1.323                        | 1.433        | 1.627         |
| q-PCR (run time app. 60 min)                                                                                                                            | DKK | 1.764                        | 1.874        | 2.068         |
| MiSeq                                                                                                                                                   | DKK | 68                           | 178          | 372           |
| NovaSeq                                                                                                                                                 | DKK | 278                          | 388          | 577           |
| Sample prep any kind of instrument for sequencing                                                                                                       | DKK | 278                          | 388          | 582           |
| Startup meeting Sequence                                                                                                                                | DKK | 0                            | 0            | 0             |
| <b>NGS instruments (service booking)</b>                                                                                                                |     |                              |              |               |
| Sequencing run start up                                                                                                                                 | DKK |                              | 2625         |               |
| SPRI beads 1ml                                                                                                                                          | DKK |                              | 441          |               |
| <b>Bioinformatics</b>                                                                                                                                   |     |                              |              | <b>Full/h</b> |
| Bioinformatics service                                                                                                                                  | DKK |                              |              | 630           |
| Startup meeting Bioinformatics                                                                                                                          | DKK | 0                            | 0            | 0             |
| <i>When you use the MSD reader you need to make a booking to reserve the instrument and a service booking to pay for the number of plates you read.</i> |     |                              |              |               |
| <b>MSD reader</b>                                                                                                                                       |     | <b>DIY/h</b>                 | <b>Min/h</b> | <b>Full/h</b> |
| Use of instrument                                                                                                                                       | DKK | 0                            | 241          | 367           |
| Each plate that is read                                                                                                                                 | DKK |                              | 278          |               |
| On-site preparation incl PBS + plastic consumables + reading (per plate)                                                                                | DKK |                              | 651          |               |
| <b>Seahorse and Cytation</b>                                                                                                                            |     | <b>DIY/h</b>                 | <b>Min/h</b> | <b>Full/h</b> |
| Use of instrument<br>(run time Seahorse app. 3h & Cytation app. 1½h)                                                                                    | DKK | 168                          | 278          | 467           |

| In Vivo Imaging System                     |     | DIY/h  | Min/h | Full/h |
|--------------------------------------------|-----|--------|-------|--------|
| Use of instrument                          | DKK | 441    | 551   | 840    |
| Services                                   |     |        |       |        |
| 1 ml Ready-to-Use aliquot (3-4 injections) | DKK | 68     |       |        |
| 1 g Potassium salt (200-300 injections)    | DKK | 4.000  |       |        |
|                                            |     |        |       |        |
| Services                                   |     | 1 hour |       |        |
| Assistance with data analyses              | DKK | 304    |       |        |